


Product details

- **Weight** : 550g
- **Power supply** : 100V 50-60Hz Power cord type Approximately 2m
- **Power consumption** : 37W
- **Frequency** : 600rpm ~ 3,600rpm
- **Frequency** : HIGH Approx. 6,000 / min LOW Approx. 5,000 / min
- **Rated time** : 30 minutes
- **Accessories** : Attachment for skin

Made in Japan



Skin Stimulation Attachment

From hard to Soft! **5** multi-head parts



⚠ Safety precautions
 ⚠ Do not use on head or neck (above bone)
 ⚠ Consult a physician before use if you are currently undergoing medical treatment or feel any physical problems.
 ⚠ Read Instruction manual carefully before use, and use only as directed.

- **Effects**
- ✓ Soothes tired muscles
 - ✓ Eases muscle stiffness
 - ✓ Relieves nerve and muscle pain
 - ✓ Improves blood circulation to help recover from fatigue

※Please note that specifications or appearance may be changed without advance notice to improve this product.

TERANISHI
 2-187 Nonami,Tenpaku-ku, Nagoya City, Aichi Prefecture 468-0045, Japan
<https://www.bullman.info>



Small yet powerful —
 the multifunctional handy massager.

5 head parts —
 choose your comfortable massage.

FighterF-86
TERANISHI

Fighter F-86



▶ *Hard & soft multi-head!*

Teranishi's mini-massager "Fighter-86" is powerful and handy-type one.

▶ *Directly pointed shiatsu!*

The head part has 5 sides, from hard to soft or even pointed- acupressure. You can choose your comfortable mode for the ideal massage.

Fighter F-86 5 effects of multi-head



● **Shiatsu surface**
Directly-pointed shiatsu massage.



● **Soft surface**
Covering wide area at once.
Soft touch massage.



● **Deep hard surface**
A massage that relieves the fatigue
of the body's deep parts.



● **Stimulation surface**
A massage that gently stimulates
the surface of the skin
to improve blood circulation.



● **Hard surface**
Covering wide area at once.
Hard touch massage.

Ways to use Fighter F-86



Neck



Shoulder



Back



Waist



Arms



Legs



Soles



Rubbing
Massage

